
Styrene C6H5CH:CH2



Performance

Measuring range	2 to 25 ppm	25 to 100 ppm	
Number of pump strokes	4(400 ml)	1 (100 ml)	
Correction factor	1	4	
Sampling time	2 min	30 sec	

Detecting limit : 0.5 ppm (4 pump strokes)

Colour change : White → Yellow Corrections for temperature & humidity : Unnecessary

Relative standard deviation: 10 % (for 2 to 5 ppm), 5 % (for 5 to 25 ppm)

Shelf life: 3 years

Reaction principle

C₆H₅CH:CH₂ + H₂S₂O₇ → Condensation polymer

Possible coexisting substances and their interferences

Substance	Concentration	Interference	Changes colour by itself to
Butadiene	≥ 5 ppm)	Dark brown
Alcohols	≥ 10 times		
Aldehydes	≥ 10 times	+ (Bleaching)	No
Esters	≥ 10 times		
Ketones	≥ 10 times	J	J

Other substance measurable with this detector tube

Substance	Correction	No. of pump strokes	Measuring range
Divinyl benzene	Factor: 0.6	3	1 to 15 ppm

Calibration gas generation

Diffusion tube method

Special note

A very low level concentration (0.2 to 4 ppm) of styrene can be measured by a Gastec special detector tube (No.124S) that is available with the Gastec Odorant Analysis System.